

Multi-Enzymatic Instrument Cleaner

(for Manual & Ultrasonic Cleaning)

Item No. 117468

Specifically developed using protease, amylase, lipase, and cellulase enzymes, and synergistically blended to produce the most powerful detergent for use on clinically-used surgical instrumentation.

Best-in-class detergent with superior soil penetration and suspension; rapidly breaks down tough-to-clean medical soils, including the multi-layers of bioburden, and prevents redeposition.

Features

- Clinically tested for use in manual cleaning and ultrasonic machines.
- Enhanced enzymatic activity combined with super detergency rapidly removes protein-rich medical soils and bioburden.
- Best in class:
 - Protease, amylase, lipase, and cellulase enzymes.
 - Corrosion inhibitors.
 - Solubilizing and suspending agents.
 - Advanced detergency and builders.
- Instantly active; Constantly active; Begins cleaning on contact.
- Low dilution rate at 1/2 oz to 1 oz per gallon (4 ml/liter – 8 ml/liter).
- Highly effective in both hard and soft water and all water temperatures.
- Safe on all critical and semi-critical medical devices; will not harm any metals, plastic, rubber, corrugated tubing, glass, or mirrors.
- Neutral pH range, non-abrasive, low-foaming, free-rinsing, non-toxic, non-corrosive, and environmentally friendly.

Directions for use

(Works in all water temperatures. 25°-30°C / 77°-86°F is optimal.)

Manual Cleaning:

- 1 Add 1/2 oz to 1 oz per gallon of water (4 ml/liter – 8 ml/liter).
- 2 Submerge instruments and appropriate scopes to be cleaned. For scopes: suction or flush through channels per scope manufacturer's instructions.
- 3 Soak for 2 minutes (up to 5 for heavily soiled instruments) to remove organic soils.
- 4 Rinse thoroughly with distilled water.
- 5 Discard the used solution.

Ultrasonic Cleaning Devices:

- 1 Set dilution rate to 1/2 oz to 1 oz per gallon (4 ml/ liter – 8 ml/liter).
- 2 Start the ultrasonic cleaner.
- 3 Penetrates instantly with superior results within 5 to 10 minutes.
- 4 Discard used solution.

Cautions:

- Avoid contact with skin and eyes. Do not swallow.
- Keep away from children.
- Wear personal protective equipment.
- Discard solution to an approved waste disposal plant.

For veterinary use only.

See SDS for any regulatory information and GHS safety directions.

To place an order, please contact your sales representative or visit mwiah.com.

Manufactured for Securos Surgical by
The Ruhof Corporation
393 Sagamore Avenue, Mineola, NY 11501
1-800-537-8463 • (516) 294-5888
Fax (516) 248-6456
www.ruhof.com • ruhof@ruhof.com
MADE IN USA